

PIPETTE FILLER - "maxi"

- A simple and easily controllable pipette filler for volumetric and graduated pipettes from 1 ml to 100 ml.
- Offers both left and right handers to operate all normal bulb and graduated pipettes easily and comfortably.
- Safe and easy handling enables even inexperienced users to adjust the meniscus precisely.
- The optimised valve system allows easy draw of liquid without excess pressure.
- Very easy to use: first squeeze the silicone bulb, move the pipetting knob up-ward for filling, down-ward for discharging. For blow-out just push the discharge bellow.
- The discharge bellow is used for blow-out (emptying) of the pipette.
- Silicone stem adaptor offers an easy insert of the pipette stem and hold the pipette in vertical position.
- A replaceable hydrophobic membrane filter in the intake cone secures cross contaminations and prevents the impurities to penetrate into the interior.

catalogue number	for pipette volume	pack quantity
011.02.003	1 ml -100 ml	1 piece



PIPETTE FILLER - "bulb"

- Bulb pipette filler features a white polypropylene attachment containing a silicone conical adapter that gently seals around pipettes without twisting.
- First squeeze the bulb, press pipette lightly into soft adapter, then relax the bulb to draw liquid.
- Pipette does not contact hard plastic which reduces the risk of cracking and chipping. Supplied with handy elastic band for hanging on reagent bottle or wrist.

catalogue number	for pipette volume	quantity in pack
011.04.001	15 ml	1 piece



PIPETTE FILLER - "bulb"

- Manufactured from pure rubber. Simple but effective tools to be used for practical filling of pipettes.
- Due to reinforced ribs they can be used with large pipettes.
- Contoured shape and heavy wall enable increased grip for easier usage.

catalogue number	suction volume	pack quantity
011.04.003	max 3 ml	10 pieces
011.04.005	max 5 ml	10 pieces
011.04.010	max 10 ml	10 pieces