



PEHANON STRIPS

- Very special pH test strip which unifies the pH indicator and the reference colors on one strip.
- pH values can be determined without a separate color chart.
- Any sample color has the same effect on both the reference colors and reactive pad.
- An invisible hydrophobic barrier just above the top color field prevents capillary action of the test solution beyond.
- Strip has 6 mm width and 85 mm length.

catalogue number	pH - range	pack quantity
MN.90415	3.8 - 5.5 pH	200 pieces



PEROXIDE STRIPS

- Peroxide strips meet the most important requirements expected from a modern quick test.
- It is a ready-to-use kit. Peroxide strips are calibrated and contains all necessary equipment and reagents.
- The user will surely receive accurate readings since the graduations are; 0 - 1 - 3 - 10 - 30 - 100.
- The color charts are adjusted using certified standard solutions which are directly tracable to primary NIST standards.
- In order to get the maximum use from the product, close the tube immediately after removing the strip.

catalogue number	strip width	strip length	pack quantity
MN.91312	6 mm	95 mm	100 pieces



AQUADUR STRIPS

- Best strip for the determination of water hardness.
- Clear color changes from green to red ensures an accurate reading. Due to the clear design of the green/red indication, a color chart is not necessary.
- The hardness of the water depends on its calcium and magnesium content.
- The total assay of these salts determines the hardness of the water.
 - < 50 ppm : very soft water
 - 50 - 120 ppm : soft water
 - 120 - 240 ppm : medium water
 - 240 - 360 ppm : hard water
 - > 360 ppm : very hard water

catalogue number	strip width	strip length	pack quantity
MN.91201	6 mm	95 mm	100 pieces

INDICATOR ROLLS - "for hot air indication"

- Perfect indicator for hot air and steam sterilisation.
- Changes color of stripe once reaches sterilisation temperature.
- Self adhesive back ensures easy and safe adhesion to almost all surfaces and textile products.

catalogue number	roll width	roll length	pack quantity
101.03.001	19 mm	50 m	1 piece

