



## DRAINING RACK

- A professional draining rack with interchangeable pegs made out of ABS.
- These push in pegs are easily removable for cleaning or to hold different shape of glassware. The peg sockets have closed ends which eliminate leakage.
  - Large size pegs are in triangular shape which can hold flasks wider than 15 mm diameter.
- 2 type of pegs:
  - Small size pegs are in round shape which can hold flasks wider than 6 mm diameter.
- The drip channel collects waste water which drains through existing drain tubes.
- Comes complete with mounting kit, draining hose and fixing template.

catalogue number	rack width	rack height	peg quantity	pack quantity
071.01.002	500 mm	610 mm	72	1 piece



NEW

## PIPETTE WASHING SET

- Washing system is ideal for gentle cleaning of pipettes based on pipetting principle. Solves the problem of washing pipettes quickly and safely.
- The polyethylene rinser fills via the spray nozzle situated on the top rim. This eliminates backsiphoning.
- Contaminated waters will syphon out totally, once the syphon point at the top of the rinser has been reached.
- Supplied as a complete system consists of a washing unit and a pipette basket.

catalogue number	quantity in pack
071.10.001	1 piece



## PIPETTE WASHING SET

- Washing system is ideal for gentle cleaning of pipettes based on pipetting principle. Solves the problem of washing pipettes quickly and safely.
- The polyethylene rinser fills via the spray nozzle situated on the top rim. This eliminates backsiphoning.
- Contaminated waters will syphon out totally, once the syphon point at the top of the rinser has been reached.
- Consists of a washing unit and cleaning vessel and a pipette basket.

catalogue number	quantity in pack
071.10.002	1 piece