


## CENTRIFUGE TUBES - "glass"

- Manufactured from neutral tubes with a wall thickness between 1,0 - 1,1 mm.
- Nominal capacity is approximately 15ml and can stand up to 3.500 RCF.
- These graduated tubes have white inscriptions on them.

| tube centrifuge | tube dimensions | wall thickness | specification | pack quantity |
|-----------------|-----------------|----------------|---------------|---------------|
| 078.01.001      | 16/17x100 mm    | 1,0/1,1 mm     | non-graduated | 100 pieces    |
| 078.01.002      | 16/17x160 mm    | 1,0/1,1 mm     | graduated     | 10 pieces     |



## TUBE CLAMP - "wood"

- Ideal for handling hot tubes or to heat test tubes up to 25 mm diameter.
- Manufactured from hardened wood with steel spring. With a long arm for better handling and to avoid close contact between fingers and the flame.
- High resistance to heat up to 70°C.

| catalogue number | total length | pack quantity |
|------------------|--------------|---------------|
| 080.04.001       | 180 mm       | 1 piece       |



## TUBE CLAMP - "st. steel"

- Ideal for handling hot tubes or to heat test tubes up to 30 mm diameter.
- Manufactured from stainless steel with a finger grip. With a long arm for better handling and to avoid close contact between fingers and the flame.
- High resistance to heat up to 100°C.

| catalogue number | total length | pack quantity |
|------------------|--------------|---------------|
| 080.04.002       | 150 mm       | 1 piece       |



## TUBE STOPPERS - "flange"

- Very effective due to their "flange" shape and grooved body.
- Available on request in different colors (red, green, yellow and blue).
- For details please refer to catalogue page 110.

