




PETRI DISHES - "glass"

- Manufactured from heat resistant glass by using counter sinking process.
- The dish and the lid are free of bubbles and the edges are fire polished.

catalogue number	diameter	height	pack quantity
081.01.060	60 mm	15 mm	18 pieces
081.01.080	80 mm	15 mm	18 pieces
081.01.100	100 mm	20 mm	18 pieces
081.01.120	120 mm	20 mm	18 pieces
081.01.150	150 mm	25 mm	18 pieces

PETRI DISHES - "polystyrene" - "sterile"



➤ Manufactured from optically clear, transparent polystyrene for multiple use in microbiology and bacteriology field.

➤ All petri dishes are manufactured under germ-free, sterile conditions of class 10.000 sterile area with extreme sterility control.

➤ Available 2 types: ➤ sterile - clean room production
➤ sterile - gamma irradiated

➤ With vented lids for consistent and effective practice. Ideal body form is suitable for automatic filling machines.

➤ 60 mm petri dishes are supplied in groups of 10 units packed in heat sealed bags.

➤ Other sizes are supplied in groups of 20 units per pack.

clean room petri dishes	diameter	height	section	pack quantity
081.02.061	60 mm	15 mm	1 section	500 pieces
081.02.091	90 mm	17 mm	1 section	500 pieces
081.02.092	90 mm	17 mm	2 section	500 pieces
081.02.093	90 mm	17 mm	3 section	500 pieces
081.02.094	90 mm	17 mm	4 section	500 pieces
081.02.121	120 mm	17 mm	1 section	500 pieces

gamma irradiated petri dishes	diameter	height	section	pack quantity
081.02.161	60 mm	15 mm	1 section	500 pieces
081.02.191	90 mm	17 mm	1 section	360 pieces
081.02.192	90 mm	17 mm	2 section	360 pieces
081.02.193	90 mm	17 mm	3 section	360 pieces
081.02.194	90 mm	17 mm	4 section	360 pieces
081.02.221	120 mm	17 mm	1 section	400 pieces