

BEAKERS - "glass" & "polypropylene"



"glass" - "low form"

- Manufactured from borosilicate glass 3.3. According to ISO 3819 and DIN 12331.
- Have thick walls and heavy bases.
- Reinforced rims and spouts to prevent breakage during washing and transportation.
- Inscriptions are in white enamel.


| catalogue number | volume | pack quantity |
|------------------|----------|---------------|
| 025.01.005 | 5 ml | 10 pieces |
| 025.01.010 | 10 ml | 10 pieces |
| 025.01.025 | 25 ml | 10 pieces |
| 025.01.050 | 50 ml | 10 pieces |
| 025.01.100 | 100 ml | 10 pieces |
| 025.01.150 | 150 ml | 10 pieces |
| 025.01.250 | 250 ml | 10 pieces |
| 025.01.400 | 400 ml | 10 pieces |
| 025.01.600 | 600 ml | 10 pieces |
| 025.01.800 | 800 ml | 10 pieces |
| 025.01.901 | 1000 ml | 10 pieces |
| 025.01.902 | 2000 ml | 6 pieces |
| 025.01.905 | 5000 ml | 1 piece |
| 025.01.910 | 10000 ml | 1 piece |



"P.P" - "embossed scale"

- Manufactured from ultra clear virgin polypropylene with superior chemical resistance according to ISO 7056.
- They have specially designed functional spouts in order to minimize possible spills.
- With embossed volume graduations and inscriptions.

| embossed scale | volume | pack quantity |
|----------------|---------|---------------|
| 026.01.025 | 25 ml | 1 piece |
| 026.01.050 | 50 ml | 1 piece |
| 026.01.100 | 100 ml | 1 piece |
| 026.01.250 | 250 ml | 1 piece |
| 026.01.500 | 500 ml | 1 piece |
| 026.01.901 | 1000 ml | 1 piece |
| 026.01.902 | 2000 ml | 1 piece |



"P.P" - "blue printed scale"

- Manufactured from ultra clear virgin polypropylene with superior chemical resistance according to ISO 7056.
- They have specially designed functional spouts in order to minimize possible spills.
- With blue printed volume graduations and inscriptions.

| blue printed scale | volume | pack quantity |
|--------------------|---------|---------------|
| 026.03.025 | 25 ml | 1 piece |
| 026.03.050 | 50 ml | 1 piece |
| 026.03.100 | 100 ml | 1 piece |
| 026.03.250 | 250 ml | 1 piece |
| 026.03.500 | 500 ml | 1 piece |
| 026.03.901 | 1000 ml | 1 piece |
| 026.03.902 | 2000 ml | 1 piece |