

BEAKERS WITH HANDLE - "stainless steel"

- Manufactured from high quality AISI 304 stainless steel for improved corrosion resistance.
- Robust seamless construction for long working life. Can withstand up to 500°C.
- Suitable for pharmaceutical and laboratory applications since all surfaces are non-reacting and non-absorbing.
- All surfaces are smooth and being highly polished enabling easy cleaning. Functional spout minimize possible spills.
- Moulded graduations on the side walls offers easy volume indication. Flared base provides an additional stability.

catalogue number	volume	diameter	height	pack quantity
025.05.500	500 ml	90 mm	140 mm	1 piece
025.05.901	1000 ml	110 mm	160 mm	1 piece
025.05.902	2000 ml	150 mm	220 mm	1 piece

**BEAKER TONGS - "chromated steel"**

- Manufactured from chromated steel. Insulated jaws are covered with special coating with perfect insulation enabling to hold all types of vessels in most secure way.

catalogue number	total length	holding capacity	pack quantity
026.03.001	260 ± 10 mm	1000 ml	1 piece



- Manufactured from chromated steel. Insulated jaws are covered with special coating with perfect insulation enabling to hold all types of vessels in most secure way.

catalogue number	total length	holding capacity	pack quantity
026.03.002	325 ± 10 mm	1000 ml	1 piece