

## ERLENMEYERS - "screw cap" - "glass" &amp; "polypropylene"



## "glass"

- Manufactured from borosilicate glass 3.3.
- Approximate volumes are printed on the body of the erlenmeyer.
- Volume graduations and inscriptions are in white enamel.
- Supplied with leakproof polypropylene screw caps.

catalogue number	neck size	volume	pack quantity
027.03.100	GL 32	100 ml	1 piece
027.03.250	GL 32	250 ml	1 piece
027.03.500	GL 32	500 ml	1 piece
027.03.901	GL 32	1000 ml	1 piece



## "polypropylene"

- Manufactured from ultra clear polypropylene.
- Approximate volumes are printed on the body of the erlenmeyer.
- Volume graduations and inscriptions are in blue enamel.
- Supplied with leakproof polypropylene screw caps.

catalogue number	neck size	volume	pack quantity
028.03.050	GL 40	50 ml	1 piece
028.03.100	GL 40	100 ml	1 piece
028.03.250	GL 52	250 ml	1 piece
028.03.500	GL 52	500 ml	1 piece
028.03.901	GL 52	1000 ml	1 piece

## STABILISING RINGS



## "for erlenmeyers"

- Vinyl coated rings are made from heavy alloy will not scratch the glass surface.
- Ideal product to prevent the erlenmeyer, flasks etc. fall over or floating in the water bath.
- Fit neatly around the neck of all sizes of erlenmeyers, flasks, filter flasks etc.

catalogue number	interior diameter	exterior diameter	pack quantity
027.11.048	48 mm	85 mm	1 piece
027.11.056	56 mm	100 mm	1 piece
027.11.061	61 mm	115 mm	1 piece
027.11.074	74 mm	129 mm	1 piece



Recommendation table for stabilising rings.

catalogue number	recommended for
027.11.048	100 ml - 500 ml
027.11.056	250 ml - 1000 ml
027.11.061	500 ml - 2000 ml
027.11.074	1000 ml - 4000 ml