

BÜCHNER FUNNELS



"porcelain"

➤ Manufactured from glazed porcelain according to DIN 12905 standards.

➤ Support plate holds the filter paper, prevents closure of the holes and allows maximum filtration speed.

catalogue number	inner diameter	pack quantity
042.01.070	70 mm	1 piece
042.01.090	90 mm	1 piece
042.01.110	110 mm	1 piece



"polypropylene"

➤ Manufactured from clear, virgin polypropylene with the top part having a perforated filter plate.

➤ The top part can be stacked with the lower part to form a vacuum seal.

➤ Upper and lower parts can be separated from each other for easy and efficient cleaning.

catalogue number	inner diameter	pack quantity
042.02.070	70 mm	1 piece
042.02.090	90 mm	1 piece
042.02.110	110 mm	1 piece

PROTECTIVE COLLARS



"rubber"

➤ Adapters made from rubber with special design that keeps the funnel and filtering flask always in center.

➤ These adapters also provide all-round protection reducing impact damages.

➤ A clean seal is created between the flask rim and the funnel during usage.

catalogue number	suitable for flask volume	pack quantity
040.05.001	100 ml	1 piece
040.05.002	250 - 500 ml	1 piece
040.05.003	1000 ml	1 piece
040.05.004	2000 ml	1 piece



catalogue number	diameter	height
040.05.001	17 mm	10 mm
040.05.002	25 mm	15 mm
040.05.003	32 mm	25 mm
040.05.004	42 mm	25 mm