

FILTER PAPERS - “quantitative”

- Perfect for quantitative routine analysis which requires high purity.
- They are acid-washed and have an extremely low ash content which is around < 0.01%.
- The total percentage of alpha-cellulose is about 95%.

catalogue number	colour code	filtration speed	diameter	pack quantity
105.01.110	blue	slow	110 mm	100 pieces
105.02.110	white	medium	110 mm	100 pieces
105.03.110	black	fast	110 mm	100 pieces
105.01.125	blue	slow	125 mm	100 pieces
105.02.125	white	medium	125 mm	100 pieces
105.03.125	black	fast	125 mm	100 pieces



FILTER PAPERS - “qualitative”

- Perfect for qualitative routine analysis and for analytical gravimetric applications.
- Wet strengthened, particularly suited for applications where the residue is removed from the filter e.g. with a spatula.

catalogue number	grade	filtration speed	diameter	pack quantity
106.01.110	No:2	slow	110 mm	100 pieces
106.02.110	No:1	medium	110 mm	100 pieces
106.03.110	No:4	fast	110 mm	100 pieces
106.01.125	No:2	slow	125 mm	100 pieces
106.02.125	No:1	medium	125 mm	100 pieces
106.03.125	No:4	fast	125 mm	100 pieces



FILTER PAPERS - “sheet”

Manufactured from high quality cellulose fibres to ensure a good resistance to humidity and a high absorption capacity.

- Widely used in general laboratory filtration, labware sterilisation, pharmaceutical product preparation rough sample preparation, also for laboratory bench protection.

catalogue number	paper dimensions	paper specifications	pack quantity
107.40.040	40 x 40 cm	85 gr / m ²	200 pieces



FORCEPS - “for filter paper”

- Manufactured from stainless steel and a perfect tool to use for handling papers.

catalogue number	pack quantity
048.04.001	1 piece

