


SPATULAS - "double end"




➤ Manufactured from high quality (AISI 304) polished stainless steel.

➤ Robust seamless construction guarantees long working life and smooth finish enables easy cleaning.

➤ Specially efficient and practical to hold products from flat and sticky surfaces due to the bent flat end.

catalogue number	spatula total length	pack quantity
047.07.130	130 ± 5 mm	1 piece
047.07.150	155 ± 5 mm	1 piece
047.07.185	180 ± 5 mm	1 piece
047.07.210	210 ± 5 mm	1 piece
047.07.250	255 ± 10 mm	1 piece
047.07.300	300 ± 10 mm	1 piece

SPATULAS - "round grooved"



➤ Manufactured from high quality (AISI 304) polished stainless steel. Robust seamless construction guarantees long working life and smooth finish enables easy cleaning.

➤ Specially efficient and practical to handle powders due to the round grooved end.

➤ One end flat spatula and the other end as round grooved.

catalogue number	spatula total length	pack quantity
047.02.150	150 ± 5 mm	1 piece
047.02.180	185 ± 5 mm	1 piece
047.02.210	210 ± 5 mm	1 piece

SPATULAS - "curved"



➤ Manufactured from high quality (AISI 304) polished stainless steel.

➤ Robust seamless construction guarantees long working life and smooth finish enables easy cleaning.

catalogue number	spatula total length	pack quantity
047.03.150	150 ± 5 mm	1 piece
047.03.180	180 ± 5 mm	1 piece
047.03.210	210 ± 5 mm	1 piece