

NSZ-405 Zoom Stereo Microscope



Microteb





Tel: 021 66 97 60 37

Fax: 021 66 90 33 54

www.microteb.com

Description:

- High quality, high function and excellent cost performance.
- Best choice with ergonomic design, sharp image display and wide viewing field.
- Providing a wide viewing field of $\Phi 20\text{mm}$ with the zoom ratio of 1:4.5.
- Ideal instrument in academic, medical and industrial field.

 <p>97mm</p>	High resolution, Wide viewing field, Long working distance
	Diopter adjustable
	Zoom objective 0.7x-3x/ 1x-4.5x, Zoom ratio 1:4.3/ 1:4.5
 <p>45°</p>	Binocular head inclined at 45°, interpupillary 52-75mm

Specification:

		NSZ-405
Drawing Tube	Binocular head, inclined at 45°	●
	Binocular head, inclined at 60°	○
	Trinocular head, inclined at 45°	○
Eyepiece	Wide field eyepiece 10×/ Φ20	●
	15×/ 16	○
	20×/ 12	○
	10× Eyepiece with cross hair	○
Objective	Zoom objective	●
	0.7×~3×	○
	1×~4.5×	○
	Auxiliary objective	○
	2×	○
Zoom Ratio	1:4.3	○
	1:4.5	●
Working Distance	97mm	●
Illumination	LED Reflection illumination	●
	LED Transmission illumination	●
Photo Attachment		○
Video Attachment		○

Note: "●" In Table Is Standard Attachment. "○" Is Optional Accessories.