



Zoom Stereo Microscope NSZ-606

Zoom Stereo Microscope NSZ-606

With high quality optical components, providing sharp and striking contrast image and ensure flat image at big depth of field.

Ease the tiredness with ergonomic design.

With LED light for both incident and transmitted illumination, providing even illumination and life expectancy can reach 6000 hours.

With Complete eyepieces, objectives and accessories, it is an ideal instrument for industry assembling, inspection and teaching field.



Technical Specifications

Drawing Tube	Binocular viewing head, Inclined at 45°
	Trinocular viewing head, Inclined at 45°
Diopter Adjustment Eyepiece	Extra wide field eyepiece WF10X/Φ22
Zoom Objective	0.8X-5X
Zoom Ratio	1:6.3
Working Distance	115mm
Glass Insert	Glass Insert, Diameter 100mm
Illumination	Incident Illumination 100V-240V/LED
	Transmitted Illumination 100V-240V/LED

