



## Stereoscopic Microscope NSZ-800

Unrivaled optical performance with infinity parallel optical system.

Providing even illumination and 6000 hours life expectancy with ultra-thin LED illumination for both incident and transmitted light.

it can be used widely for scientific research, industry, biological and medical and IVF field.

## Technical Specifications

<b>Optical System</b>	Parallel optical zoom system
<b>Total Magnification</b>	8X~80X
<b>Viewing Head</b>	Binocular head, 20°inclined, Interpupillary distance 55-75mm
<b>Eyepiece</b>	WF10X/Φ22
<b>Objective</b>	Zoom Objective 0.8X~8X
<b>Working Distance</b>	78mm
<b>Focusing Range</b>	105mm
<b>Illumination</b>	Transmission & reflection LED illumination, brightness adjustable
<b>Package</b>	Carton

