

NJC160 Metallurgical Microscope



Microteb

Tel: 021 66 97 60 37

Fax: 021 66 90 33 54

www.microteb.com

Specification:

Items	Specification	NJC-160
Viewing Head	Seidentopf type trinocular head	●
Eyepiece	EW10x/22 Extra wide field eyepiece EW10x/22	●
Infinite Plan Achromatic Objectives	4x/0.10/∞/- WD 17.8mm	●
	5x/0.12/∞/- WD 15.5mm	○
	10x/0.25/∞/- WD 10.0mm	●
	20x/0.40/∞/- WD 5.10mm	●
	40x/0.60/∞/0 WD 2.90mm	●
	50x/0.75/∞/- WD 0.32mm	○
	100x/0.80/∞/0 WD 2.00mm	○
Nosepiece	Quadruple nosepiece	●
Stage	Size 300mm×268mm, Movement range 250mm×250mm	●
Focusing	Coaxial fine & coarse adjustment, Range 24mm	●
Illumination	220V or 110V/6V20W Halogen lamp, Brightness adjustable	●
Polarization Attachment	Analyzer and polarizer	○
Filter	Blue, Green, Yellow, Ground glass	●
Photo Attachment		○
Video Adapter	C Mount	○
Camera System		○
Software	Image processing software	○

Note: "●" In Table Is Standard Attachment. "○" Is Optional Accessories.