


amber burettes

class "AS"



- Manufactured from pre-vacuum calibrated borosilicate glass.
- They meet the highest quality standards according to **DIN EN ISO 385** and **DIN 12700**.
- Permanent quality control through all steps eliminate possible deviations from nominal values.
- Calibrated **to deliver (Ex)** and the tolerances are according to **DIN/ISO Class "AS"** standards.
- **Conformity ("H") certified**, with a batch number printed on the body.
- Batch certificate can be downloaded from our website.
- Supplied with either **PTFE** or **VALVE** stopcocks.
- Inscriptions and graduations are printed with **white enamel** which is specially designed for glass use.
- Amber burettes are manufactured without schellbach stripe.

burette straight with PTFE stopcock white graduations	burette straight with VALVE stopcock white graduations	volume	tolerance	subdivision	quantity in box
018.01.010	018.02.010	10 ml	± 0.020 ml	0.02 ml	1 piece
018.01.025	018.02.025	25 ml	± 0.030 ml	0.05 ml	1 piece
018.01.050	018.02.050	50 ml	± 0.050 ml	0.10 ml	1 piece