




### beakers - "low form" - "glass"

Manufactured from borosilicate glass 3.3. According to **ISO 3819** and **DIN 12331**.

They have thick walls and heavy bases. Reinforced rims and spouts to prevent breakage during washing and transportation. Inscriptions are in white enamel.

catalogue number	volume	approximate height	Diameter	pack quantity
025.01.005	5 ml	30 mm	22 mm	10 pieces
025.01.010	10 ml	35 mm	26 mm	10 pieces
025.01.025	25 ml	50 mm	34 mm	10 pieces
025.01.050	50 ml	60 mm	42 mm	10 pieces
025.01.100	100 ml	70 mm	50 mm	10 pieces
025.01.150	150 ml	80 mm	60 mm	10 pieces
025.01.250	250 ml	95 mm	70 mm	10 pieces
025.01.400	400 ml	110 mm	80 mm	10 pieces
025.01.600	600 ml	125 mm	90 mm	10 pieces
025.01.800	800 ml	135 mm	100 mm	10 pieces
025.01.901	1000 ml	145 mm	105 mm	10 pieces
025.01.902	2000 ml	185 mm	130 mm	6 pieces
025.01.905	5000 ml	270 mm	170 mm	1 piece
025.01.910	10000 ml	325 mm	228 mm	1 piece



### erlenmeyers - "narrow neck" - "glass"

Manufactured from borosilicate glass 3.3. According to **ISO 1773** and **DIN 12380**.

Reinforced rims to prevent breaking during washing and transportation. Inscriptions are in white enamel.

catalogue number	volume	height	bottom diameter	pack quantity
027.01.025	25 ml	75 mm	42 mm	10 pieces
027.01.050	50 ml	85 mm	51 mm	10 pieces
027.01.100	100 ml	110 mm	65 mm	10 pieces
027.01.250	250 ml	145 mm	85 mm	10 pieces
027.01.300	300 ml	160 mm	90 mm	10 pieces
027.01.500	500 ml	175 mm	105 mm	10 pieces
027.01.901	1000 ml	220 mm	130 mm	10 pieces
027.01.902	2000 ml	285 mm	165 mm	6 pieces
027.01.903	3000 ml	310 mm	187 mm	1 piece
027.01.905	5000 ml	365 mm	220 mm	1 piece



### flasks - "flat & round bottom" - "glass"

Manufactured from borosilicate glass 3.3. According to **ISO 1773** and **DIN 12348**.

NS necks are in compliance with **DIN 12348** standards. Reinforced rims are used to prevent breakage during washing and transportation.

flat bottom flasks	round bottom flasks	volume	approx. diameter	pack quantity
029.01.050	029.02.050	50 ml	51 mm	10 pieces
029.01.100	029.02.100	100 ml	64 mm	10 pieces
029.01.250	029.02.250	250 ml	85 mm	10 pieces
029.01.500	029.02.500	500 ml	105 mm	10 pieces
029.01.901	029.02.901	1000 ml	131 mm	10 pieces
029.01.902	029.02.902	2000 ml	166 mm	6 pieces
029.01.904	029.02.904	4000 ml	207 mm	1 piece
029.01.906	029.02.906	6000 ml	236 mm	1 piece
029.01.910	029.02.910	10000 ml	279 mm	1 piece