

flasks - "ground neck"




Made from borosilicate glass 3.3. They meet the highest quality standards according to **DIN 12348**.

NS necks are made in compliance with **DIN 12242** standards.

Inscriptions are in white enamel which is specially designed for glass.

flasks (NS) flat bottom	flasks (NS) round bottom	volume	neck NS size	pack quantity
030.01.050	030.02.050	50 ml	NS 14/23	1 piece
030.01.051	030.02.051	50 ml	NS 19/26	1 piece
030.01.052	030.02.052	50 ml	NS 24/29	1 piece
030.01.053	030.02.053	50 ml	NS 29/32	1 piece
030.01.100	030.02.100	100 ml	NS 14/23	1 piece
030.01.101	030.02.101	100 ml	NS 19/26	1 piece
030.01.102	030.02.102	100 ml	NS 24/29	1 piece
030.01.103	030.02.103	100 ml	NS 29/32	1 piece
030.01.250	030.02.250	250 ml	NS 19/26	1 piece
030.01.251	030.02.251	250 ml	NS 24/29	1 piece
030.01.252	030.02.252	250 ml	NS 29/32	1 piece
030.01.500	030.02.500	500 ml	NS 24/29	1 piece
030.01.501	030.02.501	500 ml	NS 29/32	1 piece
030.01.900	030.02.900	1000 ml	NS 24/29	1 piece
030.01.901	030.02.901	1000 ml	NS 29/32	1 piece
030.01.902	030.02.902	2000 ml	NS 29/32	1 piece
030.01.904	030.02.904	4000 ml	NS 29/32	1 piece

erlenmeyers - "ground neck"



Made from borosilicate glass 3.3. They meet the highest quality standards according to **DIN 12387**.

NS necks are made in compliance with **DIN 12242** standards.

Inscriptions are in white enamel which is specially designed for glass.

catalogue number	volume	approx. height	neck NS size	pack quantity
028.01.050	50 ml	85 mm	NS 14/23	1 piece
028.01.051	50 ml	85 mm	NS 19/26	1 piece
028.01.052	50 ml	85 mm	NS 24/29	1 piece
028.01.053	50 ml	85 mm	NS 29/32	1 piece
028.01.100	100 ml	100 mm	NS 14/23	1 piece
028.01.101	100 ml	100 mm	NS 19/26	1 piece
028.01.102	100 ml	100 mm	NS 24/29	1 piece
028.01.103	100 ml	100 mm	NS 29/32	1 piece
028.01.250	250 ml	140 mm	NS 19/26	1 piece
028.01.251	250 ml	140 mm	NS 24/29	1 piece
028.01.252	250 ml	140 mm	NS 29/32	1 piece
028.01.300	300 ml	156 mm	NS 29/32	1 piece
028.01.500	500 ml	175 mm	NS 19/26	1 piece
028.01.501	500 ml	175 mm	NS 24/29	1 piece
028.01.502	500 ml	175 mm	NS 29/32	1 piece
028.01.901	1000 ml	220 mm	NS 29/32	1 piece
028.01.902	2000 ml	270 mm	NS 29/32	1 piece

erlenmeyers - "screw cap" - "P.P"



Manufactured from ultra clear polypropylene.

Inscriptions are in blue enamel which is specially designed for plastic.

Approximate volumes are printed on the body of the erlenmeyer.

Supplied with leakproof polypropylene screw caps.

catalogue number	neck size	volume	pack quantity
028.03.050	GL 40	50 ml	1 piece
028.03.100	GL 40	100 ml	1 piece
028.03.250	GL 52	250 ml	1 piece
028.03.500	GL 52	500 ml	1 piece
028.03.901	GL 52	1000 ml	1 piece

erlenmeyers - "screw cap" - "glass"



Manufactured from borosilicate glass 3.3.

Inscriptions are in white enamel which is specially designed for glass.

Approximate volumes are printed on the body of the erlenmeyer.

Supplied with leakproof polypropylene screw caps.

catalogue number	neck size	volume	pack quantity
027.03.100	GL 32	100 ml	1 piece
027.03.250	GL 32	250 ml	1 piece
027.03.500	GL 32	500 ml	1 piece
027.03.901	GL 32	1000 ml	1 piece