



narrow neck bottles glass

- ▶ Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- ▶ Autoclavable at 121°C and 1.1 bar for 30 minutes.
- ▶ The angled finishing allows the liquid to flow without any restrictions.
- ▶ Cap has a GL 45 thread and pouring ring which prevents drop splash.
- ▶ Amber bottles prevent from UV lights to penetrate inside the bottle.
- ▶ Inscriptions are in white enamel which is specially designed for glass.

clear bottles	amber bottles	volume	clear in box	amber in box
061.01.100	061.02.100	100 ml	10 pieces	1 piece
061.01.250	061.02.250	250 ml	10 pieces	1 piece
061.01.500	061.02.500	500 ml	10 pieces	1 piece
061.01.901	061.02.901	1000 ml	10 pieces	1 piece
061.01.902	061.02.902	2000 ml	6 pieces	1 piece
061.01.905		5000 ml	1 piece	



narrow neck bottles polyethylene

- ▶ Manufactured from high density polyethylene which can autoclavable and resist high pressures before rupture.
- ▶ The angled finishing allows the liquid to flow without any restrictions.
- ▶ Cap has a GL 45 thread and pouring ring which prevents drop splash.
- ▶ Amber bottles are manufactured from a special resin which prevent from UV lights to penetrate inside the bottle.
- ▶ Inscriptions and graduations are in embossed form.

clear bottles	amber bottles	volume	quantity in box
061.05.100	061.06.100	100 ml	10 pieces
061.05.250	061.06.250	250 ml	10 pieces
061.05.500	061.06.500	500 ml	10 pieces
061.05.901	061.06.901	1000 ml	10 pieces
061.05.902	061.06.902	2000 ml	6 pieces