

electrical  
**pipette pump**  
 adjustable nozzle



- Adjustable nozzle pipette pump combines a powerful motor with a lightweight, ergonomic design; a reflection of our commitment to the comfort of the end user.
- Features such as textured grip, hand-neutral design for right- and left-handed users, concave triggers, and convenient controls ensure the filler to perform many hours of stress-free filling and dispensing.
- Filling and dispensing speed is easily set with a single switch. Force liquid out of your pipette by setting the unit to low- or high-speed blow-out mode, or flip to gravity mode for use with "to contain" pipettes and dispense liquid without the motor.
- Filler works with glass and plastic pipettes of various volumes from 1 to 100 ml. The filter in the polypropylene nozzle prevents contamination.
- A red battery indicator warns when the batteries need to be recharged. Filler can be recharged while in use, or while resting in either the dual-position desktop or magnetic wall-mountable holder (both included).
- Pipette pump contains two additional 0.45 urn filters.

**Catalogue No: 010.01.002**

**Position One - 0°**



*The standard 0° angle (vertical position) works well for basic benchtop uses.*

**Position Two - 7 1/2°**



*The smaller-angled position works well for tight benchtops with all kind bottles and flasks.*

**Position Three - 15°**



*The 15° angled position is ideal for restricted spaces such as fume hoods and clean cabinets.*



**contains a stand that allows recharging and a magnetic holder for metal surface or wall mounting.**