



glass bottles with side arm

- ▶ Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- ▶ They are ideal products to be used with burettes and dispensers since the side arm allows to fill the bottle without removing the burette or the cap.
- ▶ These bottles have GL45 thread screw cap and a P.P pouring ring.
- ▶ Amber bottles prevent from UV lights to penetrate inside the bottle.
- ▶ Inscriptions are in white enamel which is specially designed for glass.

clear glass bottles	amber glass bottles	volume	quantity in pack
007.02.001	007.02.002	1000 ml	1 piece
007.02.003	007.02.004	2000 ml	1 piece

polyethylene bottles with side arm



- ▶ Manufactured from high density polyethylene which can resist high pressures before rupture.
- ▶ They are ideal products to be used with burettes and dispensers since the side arm allows to fill the bottle without removing the burette or the cap.
- ▶ These bottles have GL45 thread screw cap and a P.P pouring ring.

clear bottles	amber bottles	volume	quantity in pack
007.02.005	007.02.006	1500 ml	1 piece