



## flat knife spatulas

- Manufactured from excellent grade polished stainless steel such as one end wooden handle and the other end as flat knife.
- Specially efficient and practical to general use due to the flexible flat knife end.

spatula flat knife	spatula blade length	quantity in pack
047.10.075	75 mm	1 piece
047.10.100	100 mm	1 piece
047.10.130	130 mm	1 piece
047.10.150	150 mm	1 piece
047.10.165	165 mm	1 piece
047.10.200	200 mm	1 piece
047.10.250	250 mm	1 piece
047.10.300	300 mm	1 piece



P.P scoops are also suitable to be used as weighing scoops

## scoops polypropylene

- Manufactured from polypropylene. Volume marked.
- Scoops have flat bottom provides to stand horizontally.

scoops polypropylene	scoop volume	scoop length	quantity in pack
037.17.025	25 ml	135 mm	1 piece
037.17.050	50 ml	160 mm	1 piece
037.17.100	100 ml	200 mm	1 piece
037.17.250	250 ml	260 mm	1 piece
037.17.500	500 ml	315 mm	1 piece
037.17.901	1000 ml	385 mm	1 piece

## scoops stainless steel

- Manufactured from polished stainless steel without any edges and crevices for pharmaceutical use.



scoops stainless steel	scoop diameter	scoop length	quantity in pack
037.08.050	50 mm	100 mm	1 piece
037.08.075	75 mm	145 mm	1 piece
037.08.100	100 mm	190 mm	1 piece
037.08.120	120 mm	280 mm	1 piece