



tubings - "rubber" - "silicone"

Rubber tubings are manufactured from natural elastic rubber. With temperature stability -25°C to $+70^{\circ}\text{C}$.

Silicone tubings are produced biologically inert silicone. With temperature stability -60°C to $+170^{\circ}\text{C}$.

rubber tubings	silicone tubings	inner diameter	outer diameter	thickness	pack quantity
052.01.001	052.02.001	6 mm	12 mm	3 mm	1 meter
052.01.002	052.02.002	8 mm	14 mm	3 mm	1 meter
052.01.003		6 mm	15 mm	5 mm	1 meter
052.01.004		8 mm	17 mm	5 mm	1 meter



hexagonal stoppers - "hallow glass"

Manufactured from borosilicate glass 3.3, according to quality standards **DIN 12252**.

Lower closed end of stopper features a "drip-tip" to ensure adhered liquid is returned to the flask. Hexagonal flat head prevents the stopper from rolling and allows stopper to stay stable.

clear stoppers	amber stoppers	NS size	pack quantity
051.03.001	051.04.001	NS 10/19	1 piece
051.03.002	051.04.002	NS 12/21	1 piece
051.03.003	051.04.003	NS 14/23	1 piece
051.03.004	051.04.004	NS 19/26	1 piece
051.03.005	051.04.005	NS 24/29	1 piece
051.03.006	051.04.006	NS 29/32	1 piece
051.03.007	051.04.007	NS 34/35	1 piece



stopper - "rubber" - "silicone"

Rubber stoppers are manufactured from natural elastic rubber.

The temperature stability for **rubber stoppers** is between -25°C to $+70^{\circ}\text{C}$.

Silicone stoppers are produced from biologically inert silicone.

The temperature stability for **silicone stoppers** is between -60°C to $+170^{\circ}\text{C}$.

rubber stoppers	silicone stoppers	bottom diameter	top diameter	height	pack quantity
051.01.001	051.02.001	5 mm	10 mm	25 mm	10 pieces
051.01.002	051.02.002	8 mm	13 mm	25 mm	10 pieces
051.01.003	051.02.003	11 mm	16 mm	25 mm	10 pieces
051.01.004	051.02.004	14 mm	20 mm	30 mm	10 pieces
051.01.005	051.02.005	18 mm	24 mm	30 mm	10 pieces
051.01.006	051.02.006	21 mm	28 mm	35 mm	10 pieces
051.01.007	051.02.007	26 mm	34 mm	40 mm	10 pieces
051.01.008	051.02.008	32 mm	42 mm	45 mm	10 pieces
051.01.009	051.02.009	40 mm	51 mm	50 mm	10 pieces
051.01.010	051.02.010	48 mm	60 mm	55 mm	10 pieces
051.01.011	051.02.011	55 mm	68 mm	60 mm	10 pieces

