



thermometers - "solid stem"

Solid Stem Glass Thermometers are calibrated for partial immersion depth of 76 mm.

Thermometer has red alcohol column and the body is reflecting white colour which provides easier reading.

catalogue number	operation range °C	subdivision	pack quantity
059.02.001	-20/+ 50	0.5 °C	10 pieces
059.02.002	-20/+110	0.5 °C	10 pieces
059.02.003	-10/+150	1 °C	10 pieces
059.02.004	-10/+200	1 °C	10 pieces



thermometer clamp - "with bosshead"


Thermometer Clamp & Holder with boss head has an opening of 0-30 mm.

With swivel mechanism which enables the object to swivel 360° and is lockable at any desired angle.

Steel spring arm holds the thermometer firmly and keeps the thermometer at the center.

catalogue number	pack quantity
060.04.001	1 piece

thermometers - "maximum & minimum"



Thermometer with a push button for automatic resetting.

Has a weather resistant polypropylene housing with inclined shelter.

The temperature measuring range is between -30°C to +50°C, with subdivisions of 1°C.

catalogue number	pack quantity
059.03.001	1 piece

thermometers - "insertion"



Ideal for temperature measurements in liquids, powders and food products. LCD display size is 18x8 mm, probe length is 127 mm.

Operational temperature range in degree centigrade is between -50°C to +300°C, with subdivisions of 0.1°C.

Operational temperature range in degree fahrenheit is between -58 °F to +572°F with subdivisions of 0.2°F.

catalogue number	pack quantity
060.01.001	1 piece