

crucible lids

Manufactured from glazed porcelain.
Can resist temperatures up to 1000°C.

catalogue number	internal diameter	pack quantity
038.15.030	30 mm	6 pieces
038.15.035	35 mm	6 pieces
038.15.040	40 mm	6 pieces
038.15.045	45 mm	6 pieces
038.15.050	50 mm	6 pieces
038.15.060	60 mm	6 pieces



crucibles - "porcelain" - "1000°C"

Manufactured from glazed porcelain both inside and outside according to **DIN 12904** standard. Crucibles can resist temperatures up to 1000°C. **Supplied without cover.** Lids should be ordered separately.

catalogue number	external diameter	height	volume	pack quantity
038.13.H30	30 mm	38 mm	15 ml	6 pieces
038.13.M35	35 mm	28 mm	13 ml	6 pieces
038.13.H35	35 mm	44 mm	24 ml	6 pieces
038.13.L40	40 mm	25 mm	15 ml	6 pieces
038.13.M40	40 mm	32 mm	20 ml	6 pieces
038.13.H40	40 mm	50 mm	35 ml	6 pieces
038.13.L45	45 mm	28 mm	25 ml	6 pieces
038.13.M45	45 mm	36 mm	30 ml	6 pieces
038.13.H45	45 mm	56 mm	50 ml	6 pieces
038.13.L50	50 mm	32 mm	29 ml	6 pieces
038.13.M50	50 mm	40 mm	39 ml	6 pieces
038.13.H50	50 mm	65 mm	75 ml	6 pieces
038.13.L60	60 mm	38 mm	45 ml	6 pieces
038.13.M60	60 mm	50 mm	85 ml	6 pieces
038.13.M70	70 mm	60 mm	140 ml	6 pieces



crucibles - "porcelain" - "1150°C"

HALDENWANGER crucibles manufactured from glazed porcelain both inside and outside according to **DIN 12904** standard. Crucibles can resist temperatures up to 1150°C. **Supplied without cover.**



catalogue number	manufacturer cat. no.	external diameter	height	volume	pack quantity
038.15.790	79C/0	30 mm	38 mm	15 ml	5 pieces
038.15.791	79C/1	35 mm	44 mm	24 ml	5 pieces
038.15.792	79C/2	40 mm	50 mm	40 ml	5 pieces



tong - "crucible"

Manufactured from stainless steel. Curved teeth of the tong provides easy handling.

catalogue number	total length	pack quantity
038.06.250	250 mm	1 piece
038.06.300	300 mm	1 piece
038.06.400	400 mm	1 piece



tong - "evaporating dish"

Made from stainless steel and the length is 300 mm.

Unique design and curved teeth of the dish tong enables easy handling of all dishes.

catalogue number	total length	pack quantity
038.05.001	300 mm	1 piece

