

## VOLUMETRIC FLASKS - "glass" - "class A" - "clear"



- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They meet the highest quality standards according to DIN EN 1042 and DIN 12664-2.
- Individually calibrated through automatic calibration process carried out by robots.
- Permanent quality control through all steps eliminate possible deviations from nominal values.
- Calibrated to contain (In) and all the tolerances are according to DIN/ISO Class "A" standards.
- Conformity ("H", "DE-M") certified, with a batch number printed on the body.
- Batch certificate can be downloaded from our website.
- Inscriptions are printed with special blue enamel.
- Clear flasks have a diamond grinded "blue effect" meniscus ring mark for longer persistence and easy observation of meniscus line.
- All volumetric flasks are supplied with P.P stoppers.

Alternative for polypropylene stoppers a glass stopper can be selected from our catalogue page no 110.



catalogue number	volume	tolerance	approximate height	neck inner diameter	neck size	pack quantity
013.01.005	5 ml	± 0.025 ml	70 mm	10.00 mm	NS 10/19	2 pieces
013.01.010	10 ml	± 0.025 ml	90 mm	10.00 mm	NS 10/19	2 pieces
013.01.020	20 ml	± 0.040 ml	110 mm	10.00 mm	NS 10/19	2 pieces
013.01.025	25 ml	± 0.040 ml	110 mm	10.00 mm	NS 10/19	2 pieces
013.01.026	25 ml	± 0.060 ml	110 mm	12.50 mm	NS 12/21	2 pieces
013.01.050	50 ml	± 0.060 ml	140 mm	12.50 mm	NS 12/21	2 pieces
013.01.051	50 ml	± 0.060 ml	140 mm	14.50 mm	NS 14/23	2 pieces
013.01.100	100 ml	± 0.100 ml	170 mm	12.50 mm	NS 12/21	2 pieces
013.01.101	100 ml	± 0.100 ml	170 mm	14.50 mm	NS 14/23	2 pieces
013.01.150	150 ml	± 0.150 ml	200 mm	14.50 mm	NS 14/23	2 pieces
013.01.200	200 ml	± 0.150 ml	210 mm	14.50 mm	NS 14/23	2 pieces
013.01.250	250 ml	± 0.150 ml	220 mm	14.50 mm	NS 14/23	2 pieces
013.01.300	300 ml	± 0.200 ml	230 mm	14.50 mm	NS 14/23	2 pieces
013.01.400	400 ml	± 0.200 ml	250 mm	18.80 mm	NS 19/26	2 pieces
013.01.500	500 ml	± 0.250 ml	260 mm	18.80 mm	NS 19/26	2 pieces
013.01.901	1000 ml	± 0.400 ml	300 mm	24.00 mm	NS 24/29	2 pieces
013.01.902	2000 ml	± 0.600 ml	370 mm	29.20 mm	NS 29/32	2 pieces
013.01.905	5000 ml	± 1.200 ml	475 mm	34.50 mm	NS 34/35	1 piece
013.01.910	10000 ml	± 2.000 ml	570 mm	45.00 mm	NS 45/40	1 piece

150 ml, 300 ml, 400 ml and 10000 ml volumetric flasks are in addition to DIN range.