

SCOOPS - "aluminium"

- Manufactured from polished aluminium without any edges and crevices for pharmaceutical use.
- Resists to most of the chemicals and suitable for autoclaves.

catalogue number	scoop total length	pack quantity
047.18.015	85 mm	1 piece
047.18.020	125 mm	1 piece
047.18.025	145 mm	1 piece
047.18.035	185 mm	1 piece
047.18.045	250 mm	1 piece
047.18.055	290 mm	1 piece
047.18.070	350 mm	1 piece

**SCOOP - "stainless steel"**

- Manufactured from polished AISI 304 stainless steel. Suitable for all laboratory purposes and especially recommended for pharmaceutical use due to free of edges and crevices.

catalogue number	scoop total length	pack quantity
047.18.002	200 ± 5 mm	1 piece

**SPOON - "stainless steel"**

- Suitable for all laboratory purposes and especially recommended for pharmaceutical use due to free of edges and crevices.

catalogue number	spoon total length	pack quantity
047.18.001	200 ± 5 mm	1 piece

**SCOOPS - "P.P" - "sterile"**

- Manufactured from high grade polypropylene.
- Individually packed and sterilized.
- Available in two different sizes and ready to use.

catalogue number	scoop total length	pack quantity
037.28.001	115 mm	1 piece
037.29.001	150 mm	1 piece

