



FORCEPS - "general use"

- Manufactured from excellent grade polished stainless steel with blunt points and with & without 1:2 tooth.
- Available 2 types:
 - blunt point without tooth.
 - blunt point with 1:2 tooth.

forceps without tooth	forceps with tooth	forceps length	pack quantity
048.01.105		105 ± 5 mm	1 piece
048.01.115	048.02.115	115 ± 5 mm	1 piece
048.01.130	048.02.130	130 ± 5 mm	1 piece
048.01.145	048.02.145	145 ± 5 mm	1 piece
048.01.160	048.02.160	160 ± 5 mm	1 piece
048.01.200		200 ± 5 mm	1 piece
048.01.250		250 ± 10 mm	1 piece
048.01.300		300 ± 10 mm	1 piece



FORCEPS - "dissecting"

- Manufactured from excellent grade polished stainless steel with fine sharp points and sprung flutes.
- Available 2 types:
 - straight body form
 - curved body form

catalogue number	forceps body form	forceps length	pack quantity
048.07.105	sharp - straight	105 ± 5 mm	1 piece
048.07.115	sharp - straight	115 ± 5 mm	1 piece
048.07.130	sharp - straight	130 ± 5 mm	1 piece
048.07.150	sharp - straight	150 ± 5 mm	1 piece
048.07.160	sharp - straight	160 ± 5 mm	1 piece
048.07.200	sharp - straight	200 ± 5 mm	1 piece
048.08.105	sharp - curved	105 ± 5 mm	1 piece
048.08.115	sharp - curved	115 ± 5 mm	1 piece
048.08.130	sharp - curved	130 ± 5 mm	1 piece



FORCEPS - "ABS" - "disposable"

- Manufactured from ABS for disposable use with corrugated handles and serrated tips.
- Fine sharp points and straight form offer many advantages.
- Available 2 types:
 - non-sterilized
 - sterilized
- Sterilized forcepsts are **sterilized by gamma radiation and individually packed.**

catalogue number	specification	forceps length	pack quantity
048.10.110	non - sterilized	110 mm	10 pieces
048.11.110	gamma - sterilized	110 mm	1 piece