

## STOPPERS - "rubber" - "with 1 hole"

- Same technical specifications as other stoppers but with one hole drilled.

catalogue number	bottom diameter	top diameter	total height	hole diameter	pack quantity
051.21.001	10,5 mm	14,5 mm	20,0 mm	3,0 mm	1 piece
051.21.002	12,5 mm	16,5 mm	20,0 mm	3,0 mm	1 piece
051.21.003	14,0 mm	18,0 mm	20,0 mm	3,0 mm	1 piece
051.21.004	17,0 mm	22,0 mm	25,0 mm	4,0 mm	1 piece
051.21.005	18,0 mm	24,0 mm	30,0 mm	4,0 mm	1 piece
051.21.006	21,0 mm	27,0 mm	30,0 mm	5,0 mm	1 piece
051.21.007	23,0 mm	29,0 mm	30,0 mm	5,0 mm	1 piece
051.21.008	26,0 mm	32,0 mm	30,0 mm	6,0 mm	1 piece
051.21.009	29,0 mm	35,0 mm	30,0 mm	6,0 mm	1 piece
051.21.010	31,0 mm	38,0 mm	35,0 mm	7,0 mm	1 piece



## STOPPERS - "cork"

- Manufactured from non-porous natural cork with excellent quality.
- With pointed conical shape, light weight and robust. Very easy to drill through.

catalogue number	bottom diameter	top diameter	total height	pack quantity
051.23.001	3,5 mm	6,5 mm	15,0 mm	10 pieces
051.23.002	5,0 mm	9,0 mm	20,0 mm	10 pieces
051.23.003	8,0 mm	12,0 mm	20,0 mm	10 pieces
051.23.004	10,5 mm	14,5 mm	20,0 mm	10 pieces
051.23.005	12,5 mm	16,5 mm	20,0 mm	10 pieces
051.23.006	14,0 mm	18,0 mm	20,0 mm	10 pieces
051.23.007	17,0 mm	22,0 mm	25,0 mm	10 pieces
051.23.008	18,0 mm	24,0 mm	30,0 mm	10 pieces
051.23.009	21,0 mm	27,0 mm	30,0 mm	10 pieces
051.23.010	23,0 mm	29,0 mm	30,0 mm	10 pieces
051.23.011	26,0 mm	32,0 mm	30,0 mm	10 pieces
051.23.012	29,0 mm	35,0 mm	30,0 mm	10 pieces
051.23.013	31,0 mm	38,0 mm	35,0 mm	10 pieces
051.23.014	36,0 mm	44,0 mm	40,0 mm	10 pieces
051.23.015	41,0 mm	49,0 mm	40,0 mm	10 pieces
051.23.016	47,0 mm	55,0 mm	40,0 mm	10 pieces
051.23.017	50,5 mm	59,5 mm	45,0 mm	10 pieces
051.23.018	56,0 mm	65,0 mm	45,0 mm	10 pieces
051.23.019	60,0 mm	70,0 mm	50,0 mm	10 pieces
051.23.020	64,5 mm	75,5 mm	55,0 mm	10 pieces
051.23.021	71,0 mm	83,0 mm	60,0 mm	10 pieces
051.23.022	79,0 mm	92,0 mm	65,0 mm	10 pieces



## STOPPER BORER SET & SHARPENER

- Stopper borer set is manufactured from nickel plated brass. Consist of 6 pieces cork borer pipe inserted to each other.

- Complete with cleaning rod having good grip and pipes of thick gauge.

catalogue number	pipe size	set of	pack quantity
051.12.001	4 - 18	6 pieces	1 piece

## STOPPER BORER SHARPENER

- Stopper borer sharpener is manufactured from nickel plated brass.

- Consist of plated metal cone with hinged cutter blade set in slot, with wooden handel.

catalogue number	pack quantity
051.13.001	1 piece

