

## STOPPERS - "hallow glass" - "hexagonal"

- Manufactured from borosilicate glass 3.3, according to the highest quality standards DIN 12252.
- Lower closed end of stopper features a "drip-tip" to ensure adhered liquid is returned to the flask.
- Hexagonal flat head prevents the stopper from rolling and allows stopper to be stood on end preventing contamination of the ground surface.

NS size	clear stoppers	pack quantity	amber stoppers	pack quantity
NS 10/19	051.03.001	10 pieces	051.04.001	1 piece
NS 12/21	051.03.002	10 pieces	051.04.002	1 piece
NS 14/23	051.03.003	10 pieces	051.04.003	1 piece
NS 19/26	051.03.004	10 pieces	051.04.004	1 piece
NS 24/29	051.03.005	10 pieces	051.04.005	1 piece
NS 29/32	051.03.006	10 pieces	051.04.006	1 piece
NS 34/35	051.03.007	10 pieces	051.04.007	1 piece



## STOPPERS - "flange" - "polyethylene"

- Manufactured from polyethylene. Very effective due to their "flange" shape and grooved body. Due to large size flanges, they prevent leakage excellently.
- Available on request in different colors (red, green, yellow and blue).

catalogue number	for tubes diameter	colour	pack quantity
051.10.010	12/13 mm	neutral	100 pieces
051.11.010	15/16 mm	neutral	100 pieces



Very effective due to their "flange" shape and grooved body. Due to large size flanges, they prevent leakage excellently.



supplied as different colours



## FOIL DISCS - "aluminium"

- Economical and disposable covers manufactured from aluminium foil for dust-free sealing, coverage and storage of various labware such as beakers, erlenmeyers, flasks and boxes.
- Does not stick to each other since they are interleaved with tissue paper.

catalogue number	diameter	foil thickness	pack quantity
058.05.080	80 mm	0,03 mm	1000 pieces / pack
058.05.100	100 mm	0,03 mm	1000 pieces / pack
058.05.125	125 mm	0,03 mm	1000 pieces / pack
058.05.150	150 mm	0,03 mm	1000 pieces / pack

