



WATCH GLASSES - "glass"

➤ Manufactured from borosilicate glass according to DIN 12341. Fire polished edges reduces chipping and breakage.

catalogue number	diameter	pack quantity
037.01.050	50 mm	10 pieces
037.01.060	60 mm	10 pieces
037.01.070	70 mm	10 pieces
037.01.080	80 mm	10 pieces
037.01.100	100 mm	10 pieces
037.01.125	125 mm	10 pieces
037.01.150	150 mm	10 pieces
037.01.200	200 mm	10 pieces



WATCH GLASSES - "polypropylene"

➤ Manufactured from light-weight, translucent, polypropylene with high resistance to chemicals.

➤ Moulded base offers more stability and stiffness.

catalogue number	diameter	pack quantity
037.22.060	60 mm	10 pieces
037.22.080	80 mm	10 pieces
037.22.100	100 mm	10 pieces
037.22.120	120 mm	10 pieces



WEIGHING BOTTLES - "glass"

➤ Manufactured from borosilicate glass 3.3 with interchangeable lids which fits the bodies perfectly.

➤ The white labelling area is suitable for writing notes.

catalogue number	diameter	height	pack quantity
037.02.254	25 mm	40 mm	1 piece
037.02.353	35 mm	30 mm	1 piece
037.02.402	40 mm	25 mm	1 piece
037.02.305	30 mm	50 mm	1 piece
037.02.503	50 mm	30 mm	1 piece
037.02.357	35 mm	70 mm	1 piece
037.02.603	60 mm	30 mm	1 piece
037.02.408	40 mm	80 mm	1 piece
037.02.803	80 mm	30 mm	1 piece
037.02.508	50 mm	80 mm	1 piece



WEIGHING BOTTLES - "polypropylene"

➤ Manufactured from polypropylene with light weight construction. Push on lids ensure an airtight seal.

catalogue number	diameter	height	pack quantity
037.31.305	30 mm	50 mm	1 piece
037.31.403	40 mm	30 mm	1 piece
037.31.407	40 mm	70 mm	1 piece
037.31.503	50 mm	30 mm	1 piece
037.31.603	60 mm	35 mm	1 piece
037.31.609	60 mm	90 mm	1 piece
037.31.701	70 mm	120 mm	1 piece