



### WEIGHING DISHES (BOATS) - "polystyrene" - "diamond form" & "square form"

- Economical containers for weighing both liquid or powder samples in a laboratory.
- Flat bottom ensures perfect stability on counters while smooth and non-porous inner surface provides accurate pour-outs with minimal sample loss.
- They are biologically inert, contamination-free, stackable and suitable for microwave moisture analysis.
- Will stands temperature between  $-10^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$ .
- Ideal for many applications such as weighing, dispensing and storing since they are resistant to diluted acids, alkalis, alcohols, and aqueous solutions.
- Can be used also as quick freeze trays for sample material, as discard trays for broken ampules and as mixing trays.
- Manufactured from antistatic polystyrene specially usable for weighing of static effected samples. They eliminate the problems caused by static built-up when weighing fine powders.



Polystyrene weighing dishes can easily bent into pouring spouts. This function enables easy and non-spill transfer of powder and liquid samples.

diamond shape non-sterile	dish volume	dish dimentionns	dish colour	non-sterile pack quantity
037.05.005	5 ml	55 x 35 mm	white	100 pieces
037.05.030	30 ml	85 x 65 mm	white	100 pieces
037.05.100	100 ml	125 x 100 mm	white	100 pieces
037.06.005	5 ml	55 x 35 mm	black	100 pieces
037.06.030	30 ml	85 x 65 mm	black	100 pieces
037.06.100	100 ml	125 x 100 mm	black	100 pieces

diamond shape sterile	dish volume	dish dimentionns	dish colour	sterile pack quantity
037.15.005	5 ml	55 x 35 mm	white	10 pieces
037.15.030	30 ml	85 x 65 mm	white	10 pieces
037.15.100	100 ml	125 x 100 mm	white	10 pieces
037.16.005	5 ml	55 x 35 mm	black	10 pieces
037.16.030	30 ml	85 x 65 mm	black	10 pieces
037.16.100	100 ml	125 x 100 mm	black	10 pieces

square shape non-sterile	dish volume	dish dimentionns	dish colour	non-sterile pack quantity
037.09.007	7 ml	44 x 44 mm	white	100 pieces
037.09.100	100 ml	80 x 80 mm	white	100 pieces
037.09.250	250 ml	140 x 140 mm	white	100 pieces
037.11.007	7 ml	44 x 44 mm	black	100 pieces
037.11.100	100 ml	80 x 80 mm	black	100 pieces
037.11.250	250 ml	140 x 140 mm	black	100 pieces

square shape sterile	dish volume	dish dimentionns	dish colour	sterile pack quantity
037.19.007	7 ml	44 x 44 mm	white	10 pieces
037.19.100	100 ml	80 x 80 mm	white	10 pieces
037.19.250	250 ml	140 x 140 mm	white	10 pieces
037.21.007	7 ml	44 x 44 mm	black	10 pieces
037.21.100	100 ml	80 x 80 mm	black	10 pieces
037.21.250	250 ml	140 x 140 mm	black	10 pieces