

WEIGHING DISHES - "aluminium"



"for moisture determination"

- Manufactured from special grade, oil free, durable, light gauge aluminium with special design for moisture determination balances.
- Aluminium weighing dishes are with slightly tapered walls for stacking.
- Can also be used for a weighing pan and for a mixing pan enabling mixing various substances.

catalogue number	diameter	height	pack quantity
037.25.100	100 mm	8 mm	50 pieces



"with handle"

- Manufactured from special grade, oil free, durable, light gauge aluminium.
- Compressed walls for extra stability. Handles ideal for labelling or holding.
- Aluminium weighing dishes are with slightly tapered walls for stacking.
- Can also be used also as drying dish and also as mixing pan.

catalogue number	diameter	height	pack quantity
037.24.045	45 mm	12 mm	100 pieces
037.24.055	55 mm	14 mm	100 pieces
037.24.065	65 mm	16 mm	100 pieces



"general purpose"

- Manufactured from special grade, oil free, durable, light gauge aluminium
- Aluminium weighing dishes are with slightly tapered walls for stacking.
- Can also be used for a weighing pan and for a mixing pan enabling mixing various substances.

catalogue number	diameter	height	pack quantity
037.23.070	70 mm	20 mm	100 pieces
037.23.085	85 mm	37 mm	100 pieces
037.23.110	110 mm	35 mm	100 pieces