

## WEIGHING SCOOPS - "glass" - "stainless steel"

*"borosilicate glass"*

➤ Glass weighing scoops are manufactured from borosilicate glass 3.3 with flat base ensures stability on the balance pans.

➤ Hollow stem enables very easy powder transfer into narrow neck flasks.

catalogue number	max. volume	pack quantity
037.07.001	3 ml	1 piece
037.07.002	6 ml	1 piece
037.07.003	10 ml	1 piece

*"stainless steel"*

➤ Stainless steel weighing scoops are manufactured from polished stainless steel free of cavities.

catalogue number	length	pack quantity
037.04.001	95 mm	1 piece
037.04.002	117 mm	1 piece
037.04.003	128 mm	1 piece

*"weighing paper"*

➤ Manufactured from impervious, and transparent glassline with smooth surface specially usable for weighing purposes.

➤ Leaves no residues after use. Supplied as booklet of 250 sheets.

catalogue number	width	length	pack quantity
037.33.001	100 mm	100 mm	1 booklet

*"weighing brush"*

➤ Fine hair weighing brush with pointed tips and wooden handle.

➤ Excellent tool for cleaning of weighing boats, weighing dishes and analytical balances.

catalogue number	pack quantity
037.32.001	1 piece