



## BAR RETRIEVERS

- Magnetic bar retrievers are an essential tool to pick up the magnetic stirring bars.
- Constructed with a powerful magnet sealed in one end and a long stem on the other end.
- Available 2 models: retriever with solid stem and retriever with flexible stem.
- Flexible magnetic bar retriever is specially helpful to reach the flask corners.

catalogue number	retriever type	stem length	pack quantity
057.02.001	solid stem	268 mm	1 piece
057.02.002	flexible stem	330 mm	1 piece



## STIRRING RODS

- Stirring rods are very useful while stirring and mixing various solutions at laboratories.
- Available 2 types: stirring rods made of glass and stirring rods made of P.P.



### “glass rods”

- Glass stirring rods are manufactured from heat resistant neutral glass.
- Stirring rods are designed with flame polished ends for preventing accidental breakage.

catalogue number	length	diameter	pack quantity
057.03.200	200 mm	5 mm	10 pieces
057.03.250	250 mm	6 mm	10 pieces
057.03.300	300 mm	6 mm	10 pieces



### “polypropylene rods”

- P.P stirring rods are manufactured from durable polypropylene.
- Unique design of the stirring rod with one end flattened like a spatula enables effective stirring with maximum efficiency.

catalogue number	length	diameter	pack quantity
057.04.250	250 mm	6 mm	10 pieces