



## BURETTE CLAMPS

- **Steel** burette clamps are made out of steel and can hold one burette at a time.
- Can be used on support rods with 8 - 14 mm diameter. Graduations and meniscus line of the burettes are clearly visible.
- **Polypropylene** burette clamps can also hold one burette at a time.
- Polypropylene clamps can be used on support rods with 10 - 12 mm diameter. Graduations and meniscus line of the burettes are clearly visible.

<i>steel clamps</i>	<i>polypropylene clamps</i>	<i>pack quantity</i>
018.05.001	018.06.001	1 piece



## BURETTE BOTTLES

- Clear & amber bottles for automatic burettes are made out of heat resistant glass.
- Total volume capacity of the bottle is 2000 ml.
- Each bottle has a NS 29/32 grinded neck that perfectly fits to the burette.

<i>clear bottles</i>	<i>amber bottles</i>	<i>pack quantity</i>
020.07.010	020.07.025	1 piece



## BURETTE BELLOWS

- Burette bellow for automatic burettes is made from best quality of natural rubber.

<i>catalogue number</i>	<i>bellow type</i>	<i>pack quantity</i>
020.08.001	standard	1 piece
020.08.002	advanced	1 piece

