



STOPWATCH - "mechanical"

- Winding stopwatch is a sturdy, chrome plated 2-scale stopwatch.
- Case diameter is 50 mm and has a weight of 110 grams.
- Consists of stainless steel, single metallic balance, nickel based alloy hairspring, and 13 jewels fixed in the structure.
- Working temperature is between -10°C to $+40^{\circ}\text{C}$.
- Operation duration is 30 minutes and resolution is 0.2 seconds.

catalogue number	dimensions	pack quantity
055.01.001	55 mm dia.	1 piece



STOPWATCH - "electronical"

- Digital stopwatch is a light weighed and easy to use stopwatch which has a water resistant plastic body.
- Offers many functions, such as;
 - Date, time and calendar
 - Stopwatch functions
 - Alarm functions
- Run time is up to 23 hour / 59 minutes / 59 seconds.
- Operation range is 59 minutes / 59 seconds with a resolution of 1/100 seconds.

catalogue number	dimensions	pack quantity
055.02.001	63 x 76 x 20 mm	1 piece



TIMER - "mechanical"

- Designed for professional use at laboratories.
- Range is 0-60 minutes with a resolution of 1 minute.
- A clear bell ring signals at the end of the preset time. Can be reset any time.

catalogue number	dimensions	pack quantity
056.01.001	60 x 60 x 34 mm	1 piece