

THERMOMETERS - "solid stem"

➤ Solid stem glass thermometers are calibrated for partial immersion depth of 76 (+/- 5 mm) mm.

➤ Thermometer has red alcohol column and the body is reflecting white colour which provides easier reading.

catalogue number	operation range °C	subdivion	pack quantity
059.02.001	-20/+ 50	0.5 °C	10 pieces
059.02.002	-20/+110	0.5 °C	10 pieces
059.02.003	-10/+150	1 °C	10 pieces
059.02.004	-10/+200	1 °C	10 pieces


THERMOMETER CLAMP - "with bosshead"

➤ Thermometer Clamp & Holder with boss is manufactured from special alloy. The screw head has an opening of 0-30 mm.

➤ Steel spring arm holds the thermometer firmly and keeps the thermometer at the center.

➤ With swivel mechanism which enables the object to swivel 360° and is lockable at any desired angle.

catalogue number	pack quantity
060.04.001	1 piece



ELECTRONIC THERMOMETER - "insertion"

➤ Perfect for temperature measurements in liquids, powders and food products.

➤ Hygienic stainless steel integral insertion probe is protected by a plastic sleeve.

➤ Temperature can be read by digital display 19 x 8 mm.

➤ Operational temperature range in degree centigrade is between -50°C and +300°C, with divisions of 0.1°C.

catalogue number	probe length	pack quantity
060.01.001	127 mm	1 piece

