



LABORATORY BRUSH - "naturel fibre / perlon"

NEW



- Excellent tools for washing and cleaning of all type of laboratory bottles, flasks, tubes, condensers, burettes, etc.
- Brush handle is manufactured from polypropylene and offers extreme easy handling and use of the brush.
- Brush material is manufactured from highest quality of fibre/perlon for long life cycle.
- Unique handle design enables hanging of the brush to a convenient place for easy reach.
- Brush stem is manufactured from corrosion proof twisted steel wire.
- Spherical sponge top part enables access to every corner of the labware and guarantees perfect cleaning results.

catalogue number	brush code	product name & description	stem length	brush diameter	brush length	pack quantity
071.02.001	a	Burette and tubings cleaning brush	845 mm	20 mm	125 mm	1 piece
071.02.002	b	Large size bottle cleaning brush	460 mm	70 mm	130 mm	1 piece
071.02.003	c	Medium size bottle cleaning brush	460 mm	50 mm	110 mm	1 piece
071.02.004	d	Small size bottle cleaning brush	380 mm	40 mm	110 mm	1 piece
071.02.005	e	Common test tubes cleaning brush	295 mm	25 mm	100 mm	1 piece
071.02.006	f	Narrow test tubes cleaning brush	235 mm	12 mm	80 mm	1 piece