



NEW

LABORATORY BRUSH - "hinged & swivel head"

- Specially recommended for curved and cornered shaped labware such as round bottom flasks.
- Uniquely designed convex swiveling brush head gets contact with all surface of the flask including difficult curves and corners and enables thorough cleaning.
- The brush head has a slightly convex shape for fast and effective cleaning even also at heavy soiled flasks.
- Strong yet flexible bristles provide long lasting trouble free life cycle. Good grip due to the solid and long wooden handle.
- The swivel brush part fits through NS 29/32 or larger or any non tapered neck with a diameter 24 mm or larger.

| catalogue number | swivel head size | handle length | pack quantity |
|------------------|------------------|---------------|---------------|
| 071.03.001 | 21x44x24 mm | 360 mm | 1 piece |

BRUSH - "for intensive cleaning"

- Specially recommended for all applications for intensive cleaning and even for dish washing.
- Brush part consists of soft nylon bristles providing long lasting trouble free life cycle.

| catalogue number | pack quantity |
|------------------|---------------|
| 071.06.001 | 1 piece |



NEW

CLEANWARE ORGANISER

- Ensures comfort at your lab with its bin at the back which is large enough sponge, wire wool or cleaning wipes.
- Stops the mess at your lab table providing more and convenient space for use.
- The soap dispenser ensures hygienic and practical consumption of liquid detergent.
- Cleanware organizer is supplied without the brushes shown at photo.

| catalogue number | pack quantity |
|------------------|---------------|
| 071.07.002 | 1 piece |



NEW