

PIPETTES - "glass" - "class AS"

- They meet the highest quality standards according to ISO 835 - 648, DIN 12687 and DIN 12699.
- Individually calibrated through automatic calibration process carried out by robots.
- Permanent quality control through all steps eliminates possible deviations from nominal values.
- Calibrated to deliver (Ex) and the tolerances are according to DIN/ISO Class "AS" standards.
- Conformity ("H", "DE-M") certified, with a batch number printed on the body.
- Batch certificate can be downloaded from our website. Inscriptions and graduations are printed with special blue enamel.

Tips are carefully formed and bevelled to ensure proper flow rates.



All pipettes have colour coded suction stems to allow fast identification and sorting.

graduated pipettes blue scale	volume	tolerance	pack quantity
021.01.001	1 ml	± 0.007 ml	10 pieces
021.01.002	2 ml	± 0.010 ml	10 pieces
021.01.005	5 ml	± 0.030 ml	10 pieces
021.01.010	10 ml	± 0.050 ml	10 pieces
021.01.020	20 ml	± 0.100 ml	10 pieces
021.01.025	25 ml	± 0.100 ml	10 pieces
021.01.050	50 ml	± 0.200 ml	10 pieces

volumetric pipettes blue scale	volume	tolerance	pack quantity
021.05.001	1 ml	± 0.008 ml	10 pieces
021.05.002	2 ml	± 0.010 ml	10 pieces
021.05.003	3 ml	± 0.010 ml	10 pieces
021.05.004	4 ml	± 0.015 ml	10 pieces
021.05.005	5 ml	± 0.015 ml	10 pieces
021.05.010	10 ml	± 0.020 ml	10 pieces
021.05.011	11 ml	± 0.020 ml	10 pieces
021.05.015	15 ml	± 0.030 ml	10 pieces
021.05.020	20 ml	± 0.030 ml	10 pieces
021.05.025	25 ml	± 0.030 ml	10 pieces
021.05.050	50 ml	± 0.050 ml	10 pieces

