



GLOVES - "nitrile"

- Manufactured from raw material free from phthalates softeners and allergic latex proteins with approximate length of 240 mm.
- CE category 1 for minimal risks for usage where both wearing comfort and product/process safety have priority. Tested in accordance with EN Norms EN420:2003.
- AQL 1,5 offers the highest level for EN374-2:2003 protection against penetration by micro organisms.
- Very good mechanical strength with excellent fit and excellent good grip through textured fingers.
- %40 stronger and %60 leaner than latex gloves.
- with standard thickness of minimum 0,06 mm.
- Available 2 types:
 - with extra thickness of minimum 0,12 mm.
- Extra thick gloves are specially designed for life sciences, biomedical research, forensic science, R&D and non-sterile drug manufacturing.
- Provides protection against wide range of chemicals including cytotoxic drugs.
- Easy to wear and particularly skin friendly through special finishing process.

standard gloves	glove size	colour code	AQL	glove length	pack quantity
080.15.001	Small	blue	1.50	240 mm	100 pieces
080.15.002	Medium	blue	1.50	240 mm	100 pieces
080.15.003	Large	blue	1.50	240 mm	100 pieces
080.15.004	X Large	blue	1.50	240 mm	100 pieces

extra thick gloves	glove size	colour code	AQL	glove length	pack quantity
080.16.001	Small	blue	1.50	240 mm	100 pieces
080.16.002	Medium	blue	1.50	240 mm	100 pieces
080.16.003	Large	blue	1.50	240 mm	100 pieces
080.16.004	X Large	blue	1.50	240 mm	100 pieces