

GLOVES - "cold resistant"

- Excellent gloves to work in cold treatment plants, with coolants and as well as for handling dry ice and similar applications.
- Offers protection against the danger when working with ultra low temperatures down to -140°C.
- Even after working with extreme temperatures they still stay flexible.

catalogue number	glove size	pack quantity
080.10.001	universal	1 pair

**GLOVES - "chemical resistant"**

- Made from polychloroprene which is resistant to many chemical substances such as acids, bases and organic oils and fats as well as bacteriological contamination.
- Anatomical shape with good gropeability with cotton flocked.
- These gloves have safety category 3 and can be used in laboratories, trade and industry.

catalogue number	glove size	glove length	pack quantity
080.12.001	Small	310 mm	1 pair
080.12.002	Medium	310 mm	1 pair
080.12.003	Large	310 mm	1 pair
080.12.004	X-Large	310 mm	1 pair

**GLOVES - "cotton"**

- Made from lightweight 100% cotton knit.
- Can be worn alone or pair and suitable for both men and women.
- Provides a comfortable and soft usage and fits either hand.

catalogue number	glove size	pack quantity
080.13.001	universal	1 pair

**GLOVES - "heat resistant"**

- Made from soft, pliant terry cotton which is machine washable.
- Provides excellent heat protection up to 250°C during autoclaving and furnacing processes.
- Very flexible and easy to wear and offers comfort in applications.

catalogue number	glove size	pack quantity
080.11.001	universal	1 pair

