

## STAINING JARS - "Hellendahl" - "Schifferdecker &amp; Papanicolau"

**"Hellendahl"**


- Hellendahl staining jar is suitable for holding 76x26 mm slides vertically.
- Hellendahl staining jar has 9 grooves which can accommodate 18 slides one after the other.
- Staining jar supplied with glass lid.

catalogue number	pack quantity
073.01.001	1 piece

**"Schifferdecker"**

- Schifferdecker staining jar is suitable for holding 76x26 mm slides horizontally.
- Schifferdecker staining jar has 10 grooves which can accommodate 20 slides one after the other.
- Staining jar supplied with glass lid.

catalogue number	pack quantity
073.03.001	1 piece

**"macro"**

- Macro staining jar is manufactured from thick moulded glass.
- Stainless steel staining rack with 20 grooves which can accommodate 20 single slides.
- Staining jar supplied with glass lid.

catalogue number		pack quantity
073.05.001	jar	1 piece
073.06.001	rack	1 piece

**NEW****"staining jars" - "PMP"**

- Made from polymethylpentene (PMP) and included lid.
- Hellendahl staining jar is suitable for holding slides vertically and Schifferdecker staining jar is suitable for holding slides horizontally.
- Remark: Xylene can cause etching to PMP after two days at room temperature.

catalogue number	model	pack quantity
073.01.011	Hellendahl	4 pieces
073.03.011	Schifferdecker	4 pieces

**NEW**