



PIPETTE STAND - "carrousel"

- Manufactured from moulded polypropylene and can hold up to 94 pipettes vertically.
- The round carrousel stand can rotate 360° around stand axis enabling very easy reach of any pipette.
- Upper disc keeps the pipettes in the correct position.
- The lower rotating disc consists of four concentric troughs with tapered sides and open bottoms which cradle the tips protectively.
- The upper rotating disc has a grid pattern concentric with the bottom disc keeping the pipettes in vertical position.
- Liquid draining of the pipette tip is retained on top of the stationary base by the outer rim.
- A low center of gravity stabilizes the 48 cm high stand on the 23 cm base diameter.

catalogue number	diameter	height	pack quantity
022.01.001	230 mm	480 mm	1 piece



PIPETTE STAND - "linear"

- Manufactured from moulded polypropylene and can hold up to 12 pipettes horizontally.
- Offers ideal storage of pipettes without any pressure on fragile pipette tips.
- Easy clean and easy dismantle for space saving storage.

catalogue number	dimensions	pack quantity
022.02.001	240 x 180 x 290 mm	1 piece



PIPETTE TRAY

- Made from moulded polypropylene. Can be used in drawers or on the bench top for storing pipettes.
- Scalloped rims can hold 7 large pipette or 16 small pipettes while the center part can hold 12 short pipettes in perpendicular direction.

catalogue number	dimensions	pack quantity
022.03.001	283 x 216 x 40 mm	1 piece