

## MICROSCOPE COVER GLASSES

- Manufactured from high optical grade, hydrolytic class no:1 borosilicate glass with an excellent chemical resistance.
- Their fully automatical production ensures clean, dust free cover glasses with accurate cut edges free of flaws and blisters.
- Thickness of the glass used for production is no:1 grade (0,13 -0,17 mm) with refractive index of n=1,52
- Available 2 types:
  - economical quality
  - superior quality

economic quality	superior quality	dimensions	quantity in pack	quantity in box
075.00.001	075.01.001	18 x18 mm	100 pieces	10 pack
075.00.002	075.01.002	20 x 20 mm	100 pieces	10 pack
075.00.003	075.01.003	22 x 22 mm	100 pieces	10 pack
075.00.004	075.01.004	24 x 24 mm	100 pieces	10 pack
075.00.005	075.01.005	24 x 32 mm	100 pieces	10 pack
075.00.006	075.01.006	24 x 50 mm	100 pieces	10 pack
075.00.007	075.01.007	24 x 60 mm	100 pieces	10 pack



## MICROSCOPE SLIDES

- Microscope slides are manufactured from high optical grade soda lime glass with hydrolytic class 3.
- Their fully automatic production ensures clean, dust and grease free slides.
- They have extremely flat surface with 1 mm thickness and 26 mm x 76 mm size.
- Available 2 types:
  - economical quality
  - superior quality



Superior quality colour coded slides have ground edges, with 45° edge angle, with 5 different colours for easy identification.



economic quality	superior quality	specification	edge angle	colour code	quantity in pack	quantity in box
075.05.001		Cut edges, Plain	90° Degree	-----	50 pieces	50 pack
075.05.002		Cut edges, Frosted one end one side	90° Degree	-----	50 pieces	50 pack
075.05.003		Ground edges, Plain	90° Degree	-----	50 pieces	50 pack
075.05.004		Ground edges, Frosted one end two side	90° Degree	-----	50 pieces	50 pack
	075.02.005	Ground edges, Colour Frosted	45° Degree	white	50 pieces	50 pack
	075.02.006	Ground edges, Colour Frosted	45° Degree	yellow	50 pieces	50 pack
	075.02.007	Ground edges, Colour Frosted	45° Degree	green	50 pieces	50 pack
	075.02.008	Ground edges, Colour Frosted	45° Degree	pink	50 pieces	50 pack
	075.02.009	Ground edges, Colour Frosted	45° Degree	blue	50 pieces	50 pack