

TEST TUBE RACKS - "P.P" - "dismantled"

- Made of "modified" polypropylene which enables the rack to be used in water baths without the risk of floating.
- Alpha numeric grid reference on top tier of the rack ensures an easy sample location.
- Racks are delivered dismantled in a bag with assembling instruction included.
- Supplied dismantled and can be assembled in less than a minute.


catalogue number	tube diameter	tube quantity	tube array	pack quantity
079.01.001	13 mm	90 pcs	6 x 15	1 piece
079.01.002	16 mm	60 pcs	5 x 12	1 piece
079.01.003	20 mm	40 pcs	4 x 10	1 piece
079.01.004	25 mm	24 pcs	3 x 8	1 piece
079.01.005	30 mm	21 pcs	3 x 7	1 piece




TEST TUBE RACKS - "P.P" - "monoblock"

- Made of "modified" polypropylene which enables the rack to be used in water baths without the risk of floating.
- Round hole tube sockets hold the tubes firmly at exact position.
- Alpha numeric grid reference on top tier of the rack ensures an easy sample location of the tubes.
- Supplied assembled and can be dismantled in less than a minute.

catalogue number	tube diameter	tube quantity	tube array	pack quantity
079.03.001	12-13 mm	50 pieces	5 x 10	1 piece
079.03.002	15 - 16 mm	50 pieces	5 x 10	1 piece
079.03.003	17 - 18 mm	50 pieces	5 x 10	1 piece
079.03.004	18 - 20 mm	45 pieces	5 x 9	1 piece
079.03.005	28 - 30 mm	26 pieces	4 x 6 / 7	1 piece



Supplied dismantled and can be assembled in less than minute.



Supplied assembled and can be dismantled in less than minute.