




TUBE RACK - "polypropylene" - "snapable"

- Manufactured from polypropylene. Designed for stable use in water baths.
- Very easy to assemble, system snaps together securely.
- Each tube row is alphanumericly matrixed for quick identification of tubes.
- Rack can hold 30 pieces of 15 ml tubes and 20 pieces of 50 ml tubes.

catalogue number	tube volume	tube quantity	tube array	pack quantity
079.04.001	15 & 50 ml	30 & 20	assorted	1 piece



Holes at the bottom plate guarantees firm hold of tubes in their position.



TUBE RACK - "polypropylene" - "4 side"

- Manufactured from polypropylene, offering a unique system of tabs and slots.
- Four sides offers to holding 4 different tube combinations:
 - Side (1); holds 8 pieces of 30 mm dia. tubes
 - Side (2); holds 20 pieces of 20 mm dia. tubes
 - Side (3); holds 24 pieces of 17 mm dia. tubes
 - Side (4); holds 40 pieces of 12 mm dia. tubes
- Racks can be linked together to form a multiple rack system.

catalogue number	tube diameter	tube volume	tube quantity	tube array	pack quantity
080.01.001	30 mm	50 ml	8	2 x 4	1 piece
	20 mm		20	4 x 5	
	17 mm	15 ml	24	3 x 8	
	12 mm		40	4 x 10	

