


RACK - "for micro tubes" - "double side"

- Manufactured from polypropylene for alternative storage of micro tubes.
- Stand has two sides and each side accommodates different size and combination of microcentrifuge tubes.
- First side can hold 96 pieces of 1,5 ml tubes, and the other side holds 48 pieces of 0,5 ml tubes and also 96 pieces of 0,2 ml tubes.
- Supplied with a clear polypropylene lid which covers the rack and protects the tubes from dust and other effects.
- Alphanumerical index on tube rows enables easy identification of tubes.

catalogue number	tube volume	tube qty.	tube array	pack quantity
080.06.002	1,5 ml	96	6 x 16	1 piece
	0,5 ml	48	6 x 8	
	0,2 ml	96	8 x 12	



First side can hold 96 pieces of 1,5 ml tubes, and the other side holds 48 pieces of 0,5 ml tubes and also 96 pieces of 0,2 ml tubes.



LABELS - "for micro tubes & cryo vials"

- Manufactured from flexible polyofiber label material with acrylic adhesive. This unique combination expands and contracts together with temperature fluctuations resulting in a strong adhesive bond.
- Chemically inert, resists most organic solvents and caustic agents. They adhere to most plastics, glass and metals without cracking, peeling or degrading.
- Pre-cut, peel-off labels which accept writing from most lab markers including all solvent resistant pens.
- Perfectly sized to fit cryo vials and cryo tubes. Withstand from -96°C to +150°C

catalogue number	tag size length x width	reccomended to be used for	roll quantity	pack quantity
099.05.033	33 x 13 mm	1,5 - 2,0 ml tubes	1000 tags	1 piece

