

MICROPIPETTE STANDS - "linear"

- Manufactured from robust and virgin polypropylene.
- Offers very easy direct access to the micropipettes. Can be used with most common all pipette brands.
- Can hold up to 6 single or multi channel pipettes in a linear form.

catalogue number	base dimensions	total height	pack quantity
006.11.006	143 x 290 mm	290 mm	1 piece



MICROPIPETTE STANDS - "adjustable"

- Manufactured from white acrylic with ergonomic angled body form and adjustable pipette holding design.
- Modern and ergonomic design enables extremely easy holding and handling of 6 single channel micropipettes.

catalogue number	base dimensions	total height	pack quantity
006.12.004	126 x 260 mm	240 mm	1 piece
006.12.006	126 x 380 mm	240 mm	1 piece



With the help of adjustable pipette holding design long and short height pipettes can be easily stored.



MICROPIPETTE STANDS - "with bin"

- Manufactured from white acrylic with clear front panel for easy identifications of bin content.
- The top shelf can hold 4 single channel micropipette.
- Self closing sliding lid guarantees dust-free storage.

catalogue number	base dimensions	total height	pack quantity
006.05.001	200 x 178 mm	173 mm	1 piece



Sliding cover with generous openings make filling and dispensing easy. This storage cabin is also ideal for dust-free storage of various types of products.

